

Lysozyme reaction buffer

Revision date: 2017/11/01

Version: 2017-01

1. PRODUCT AND COMPANY IDENTIFICATION

Product number: Lysozyme reaction buffer

Recommended use: Laboratory Chemicals

Supplier: RBC Bioscience Corp.

3F, 3F.-1,2, No. 132, Ln. 235, Bao Quao Rd., Xindian Dist.,
New Taipei City 23145, Taiwan (R.O.C.)

Emergency/Phone Number: +886-2-8912-1200

E-mail: infor@rbcbioscience.com

2. HAZARDS IDENTIFICATION

GHS Classification

Not a hazardous substance or mixture

GHS Label element

Not a hazardous substance or mixture

Potential Health Effects:

Aggravated Medical Condition: None known.

Symptoms of Overexposure: No information available.

Carcinogenicity

IARC

No ingredient of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.

ACGIH

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

OSHA

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.

NTP

No ingredient of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.

3. COMPOSITION/INFORMATION ON INGREDIENTS

Substance/Mixture: Mixture

Hazardous ingredients

No hazardous ingredients

4. FIRST AID MEASURES

General advice:

Show this material safety data sheet to the doctor in attendance..

If inhaled:

Move to fresh air. If symptoms persist, call a physician.

In case of skin contact:

Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. If symptoms persist, call a physician.

In case of eye contact:

Remove contact lenses. Protect unharmed eye. Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed:

If accidentally swallowed obtain immediate medical attention. Rinse mouth with water. Never give anything by mouth to an unconscious person.

Most important symptoms and effects, both acute and delayed:

No information available.

Note to Physician:

No information available.

5. FIRE FIGHTING MEASURE

Suitable Extinguishing media:

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Specific hazards during firefighting:

Exposure to decomposition products may be a hazard to health.

Hazardous combustion products:

No hazardous combustion products are known.

Specific extinguishing methods:

In the event of fire and/or explosion don't breathe fumes.

Special protective equipment for fire-fighters:

Wear self-contained breathing apparatus for firefighting if necessary

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures:

Use personal protective equipment. Avoid dust formation. Avoid breathing dust/fume/gas/mist/vapors/spray.

Methods and materials for containment and cleaning up:

Keep in suitable, closed containers for disposal.

7. HANDLING AND STORAGE

Advice on safe handling:

For personal protection see section 8. Smoking, eating and drinking should be prohibited in the application area.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Ingredients with workplace control parameters

Contains no substances with occupational exposure limit values.

Personal protective equipment

Hand protection remarks:

The choice of an appropriate glove does not only depend on its material but also on other quality features and is different from one producer to the other. Take note of the information given by the producer concerning permeability and break through times, and of special workplace conditions (mechanical strain, duration of contact)

Eye protection:

Safety glasses.

Skin and body protection:

Choose body protection according to the amount and concentration of the dangerous substance at the work place. Footwear protecting against chemicals.

Hygiene measures:

Keep away from food and drink. When using do not eat, drink or smoke.

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance: Liquid

Color: colorless.

Odor: characteristic.

Odor Threshold: No data available.

pH: 6-8.

Melting point/range: 0°C.

Boiling point/range: 100°C/No data available

Flash point: Not applicable.

Evaporation rate: No data available.

Burning rate: No data available.

Upper explosion limit: No data available.

Lower explosion limit: No data available.

Vapor pressure: No data available.

Relative vapor density: No data available.

Relative density: No data available.

Density: 1.0 g/cm³.

Solubility: No data available.

Water solubility: completely miscible.

Solubility in other solvents: No data available.

Partition coefficient n-octanol/water: No data available.

Auto ignition temperature: doesn't ignite..

Thermal decomposition: No data available.

Viscosity: No data available.

Viscosity, dynamic: No data available.

Viscosity, kinematic: No data available.

Explosive properties: No data available.

Oxidizing properties: No data available.

10. STABILITY AND REACTIVITY

Reactivity: No decomposition if stored and applied as directed.

Chemical stability: No decomposition if stored and applied as directed.

Possibility of hazardous reactions: Stable under recommended storage conditions.

Conditions to avoid: Heat, Flames and sparks.

Incompatible materials: No data available

Hazardous decomposition products: No decomposition if stored and applied as directed.

11. TOXICOLOGICAL INFORMATION

Acute toxicity

Product:

Acute oral toxicity: No data available

Acute inhalation toxicity: No data available

Acute dermal toxicity: No data available

Skin corrosion/irritation

Product:

May cause skin irritation susceptible persons.

Serious eye damage/eye irritation

Product:

Remarks: May irritate eyes.

Respiratory or skin sensitization

No data available

Germ Cell mutagenicity

No data available

Carcinogenicity

No data available

Reproductive toxicity

No data available

STOT-single exposure

No data available

STOT-repeated exposure

No data available

Aspiration toxicity

No data available

Further information

No data available

12. Ecological Information

Eco toxicity

Product

Toxicity to fish: No data available

Toxicity to algae: No data available

Toxicity to bacteria: No data available

Persistence and degradability

No data available

Bio accumulative potential

Product

Bio accumulation: No data available

Partition coefficient n-octanol/water: No data available

Mobility in soil

No data available

Other adverse effects

No data available

13. DISPOSAL CONSIDERATIONS

Disposal Methods

Waste from residues

Dispose of contents/container in accordance with local regulation.

Contaminated packaging:

Empty containers should be taken to an approved waste handling site for recycling or disposal containers.

14. TRANSPORT INFORMATION

UNRTDG

Not regulated as a dangerous good

IATA-DGR

Not regulated as a dangerous good

IMDG-Code

Not regulated as a dangerous good

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

No data available

Domestic regulation

49 CFR

Not regulated as a dangerous good

Special precaution for user:

Remarks: Not classified as dangerous in the meaning of transport regulations.

15. REGULATORY INFORMATION

OSHA Hazards: No OSHA Hazards

EPCRA-Emergency Planning and Community Right-to-Know

SARA 304 Extremely Hazardous Substances Reportable Quantity

This material does not contain any components with a section 304 EHS RQ.

SARA 311/312 Hazards: No SARA Hazards

SARA 302:

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313:

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

US State Regulations

California Prop 65

This product does not contain any chemicals known to the State of California to cause cancer, birth or any other reproductive defects.

16. OTHER INFORMATION

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.